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REF : 260105

INTRODUCTION

This product is manufactured from donor red blood cells with characterized blood groups and extended phenotypes. It is composed of 5 red blood cells selected for their specific phenotypic diversity across the Rhesus, Kell, Duffy, Kidd, MNSS, Lewis, P and Lutheran systems.

These "raw material" red blood cells are not ready-to-use. They are intended to be used in any techniques requiring red blood cells with known extended phenotypes.

PRESENTATION – CONTENT

Concentrated red blood cells are supplied as a kit containing 5 red blood cells. Each red blood cells unit are packed in tubes containing 1.5mL of a 60% theoretical concentration suspension in Saline-Adenine-Glucose-Mannitol (SAGM) storage solution, with the addition of a combination of antibiotics (Sulfamethoxazole and Trimethoprim) to prevent bacterial contamination.

An antigram detailing the antigenic profile of the red blood cells is provided with each kit.

CONSERVATION - VALIDITY

These samples must be stored with a cap at a temperature between +2°C and +8°C.

Their use must respect the expiry date indicated on each label, provided the product has not been subjected to any further processing.

Do not freeze.

QUALITY CONTROL

The product undergoes systematic batch release testing to ensure batch to batch consistency with quality requirements.

SYMBOLS USED

 REF	Catalogue reference		Legal manufacturer	 LOT	Batch number
	Use until		Temperature limits		

OPERATING INSTRUCTIONS – CONDITION OF USE

This product is manufactured from erythrocyte concentrate obtained from voluntary blood donors and is intended to be used as a raw material. The user is responsible for developing the application method and validating the treatment applied to the red blood cells.

This product is for *in vitro* use only.

RESULTS INTERPRETATION

Since no specific application technique is claimed for these red blood cells, the expected result is defined by the user. This outcome must take into account the extended phenotypes listed on the antigram provided with each batch.

PRECAUTION FOR USE

These reagents are prepared from whole blood collected from voluntary blood donors in accordance with ethical and technical regulations in France. They were screened for transmissible diseases in accordance with the standards in force set forth by French law. However, despite the negative result of these tests, these samples must still be considered as potentially contaminants and therefore must be handled with the usual precautions as they are biological products of human origin. Waste disposal must be carried out in accordance with current legislation.

For technical assistance or request for information, please contact us by phone:
Tel: +33 6-62-99-46-76 or e-mail: efs.reactifs.hotline@efs.sante.fr

Reference: UPR/PUL/CED/DF/FO/011
Version: 7 - 21/05/2024