**REF 232614 Hemo DTT 8**

**Principe**

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REF **232614** **Hemo DTT 8**

#### Principle

This ready-to-use product is a reducing agent from the thiol family that breaks the disulphide bonds in molecules. Treatment with 0.2 M DTT alters or even destroys certain erythrocyte antigens, the structure of which has a lot of disulphide bonds. This treatment:

- destroys antigens from the KEL, LW, SC, IN, JMH systems and the YT1 antigen
- weakens the antigens of the KN, LU, DO and AN/WJ systems.

The DTT is dissolved in Phosphate Buffer (PB) pH 8.

#### Presentation

The product consists of a range comprising 1 tube of 0.2 M DTT. The solution is packaged in 1 mL vials.

**Strictly for *in vitro* use only.**

#### Storage conditions and shelf life

This product must be stored frozen at a temperature of -20°C or below. It must not be used after the expiry date stated on the label.

Use immediately after thawing.

Beyond this point, DTT is no longer stable and non-specific and aberrant reactions may occur.

Do not refreeze the products after thawing.

If the integrity of a vial is compromised, do not use it and dispose of it in accordance with the regulations in force.

#### Quality control

The preparation of this product is subject to systematic batch controls to ensure consistent quality from batch to batch.

#### Samples

Red blood cells from a panel of test red blood cells in storage solution or red blood cells from donor tubes.

#### Conditions of use and limitations

- Check visually that there is no precipitate or substance in the DTT vial.
- Using this technique, it is possible to prepare informative red blood cells in order to provide guidance in the event of "pan-agglutination" during the identification stage of irregular antibody testing (IAT) with a negative autologous control (suspicion of alloantibodies directed against a high-frequency antigen). The preparation of these red blood cells makes it possible to focus on the public specificities of the above-mentioned systems, which will be destroyed or weakened after treatment.

**Additional equipment required:** incubator at 37°C ± 2°C, dry glass tubes, 0.15M NaCl

#### Protocol:

##### Selection of red blood cells for treatment

- Select 2 red blood cells (or 1 if necessary) from the panel of test red blood cells or donors (or patients if necessary) with group O whose RH-KEL1 phenotype is antigen-compatible with the associated antibodies (present or not eliminated).
- Ensure that the blood compatibility test (BCT) with the plasma/serum under study is positive with each selected red blood cell and perform an indirect antiglobulin test (IAT) using the tube AND filtration technique before any treatment of the red blood cells.

##### Method

- Thaw 1 tube of 0.2M DTT pH 8 to thaw at room temperature.
- Prepare 1 volume of washed red cell concentrate in 0.15 M NaCl of the selected red blood cells
- To this volume of washed red cell concentrate, add 4 volumes of DTT pH 8 0.2M at room temperature
- Incubate the mixture for 30 minutes at +37°C
- Wash 4 times in 0.15 M NaCl
- Prepare a 3-5% suspension in 0.15 M NaCl.

Quality control of the red blood cell treatment: Phenotyping the red blood cells before and after treatment in the KEL system (KEL2 or KEL4 antigen). Ensure compliance with quality control

- Perform a BCT in an IAT tube between each treated red blood cell and the plasma/serum of the patient under investigation  
The BCT must be performed in accordance with the laboratory's standard operating procedure. Obtaining positive or negative reactions will provide guidance on the specificity of the antibody in the event of "pan-agglutination".

#### Control and limitations of the technique

The use of DTT pH8 for the treatment of red blood cells to obtain diagnostic guidance is validated for the selective alteration of certain antigens, in accordance with data in the literature.

Validation of DTT treatment requires performing a KEL2 (k) phenotype (for example) before and after treatment of the selected red blood cells. Quality control is compliant if the KEL2 phenotype is positive prior to DTT treatment and negative after treatment.

#### Interpretation of the results

The results are interpreted in accordance with the table below.

BCT (IAT tube method)	RBC not treated with DTT	RBC treated with 0.2M DTT pH8	Conclusion:
Patient serum/plasma	Positive	Positive	Antibodies recognising an antigen not destroyed by DTT 0.2 M pH 8
Patient serum/plasma	Positive	Negative	Antibodies recognising an antigen destroyed by DTT 0.2 M pH 8
Patient serum/plasma	Positive	Reduced intensity	Antibodies recognising an antigen weakened by DTT 0.2 M pH 8

#### Performances

Stored at a temperature of -20°C or below and used in accordance with the technical procedures described above, the efficacy of the DTT product was verified by performing a BCT in an IAT tube between treated red blood cells and different serums from patients carrying antibodies of known specificity.

For the serums:

- Serums from patients with alloantibodies of known specificity (confirmed by identification by the CNRGS)
- Anti-RH1 from the CQI

For red blood cells treated with DTT pH 8 0.2 M

- Red blood cells from donor tubing (RBD tubing)
- Red blood cells from master solutions of red blood cells for screening panels

The BCTs carried out gave the following results:

	Specificity	BCT in IAT Tube BEFORE treatment	BCT in IAT Tube AFTER treatment
Antibodies directed against Antigens altered by 0.2 M DTT pH 8	Anti-KN Anti-AnWj Anti-Lub	Positive	Reduced intensity or even negative if Ac is weak
Antibodies directed against Antigens destroyed by 0.2 M DTT pH 8	Anti-JMH Anti-KEL2 Anti-KEL4 Anti-Yta	Positive	Negative
Antibodies directed against Antigens not altered by 0.2 M DTT pH 8	Anti-D	Positive	Positive

#### Precautions for use

The DTT in Phosphate Buffer solution is deemed less harmful than DTT before reconstitution; however, it can cause skin burns or eye damage and should therefore always be handled with the usual precautions. If in doubt, or if skin or eye irritation persists, seek medical advice. In case of contact with the eyes, rinse thoroughly with water, keeping the eyelids open. In case of contact with skin or hair, remove contaminated clothing and rinse with water. If swallowed, rinse mouth thoroughly with water (if the person is conscious). Do not induce vomiting. Hospitalize.

Waste disposal must be carried out in accordance with the legislation in force. In the event of accidental spillage, ventilate the area, avoid contact with skin or eyes, do not breathe vapours. Dilute the spilled product thoroughly with water and prevent contamination of the soil and water. Absorb the spilled product with an inert absorbent material (sand, vermiculite, etc.). Rinse the area thoroughly with water and treat the contaminated water in accordance with the regulations in force.

#### Bibliography

- Jill R.Storry and al "A Review: modification of the red blood cell membrane and its application in blood group serology. Immunohematology, vol.16,2000
- T. Tamada and T. Mazda "Enzyme and DTT treatment of adherent RBCs for antibody identification by a solid phase immunoassay system" Immunohematology, vol 18, 2002
- M. Reid The blood group Antigen. Facts Book. Second edition 2007,2008

#### Symbols used:

 REF	Catalogue reference	 LOT	Batch number		Use by		Temperature limits
	Manufacturer						

For technical assistance, information requests or complaints, please contact us at:

Tel: **06.62.99.46.76**

E-mail: [efs.reactifs.hotline@efs.sante.fr](mailto:efs.reactifs.hotline@efs.sante.fr)

UPR/PUL/FAB/DF/FI/147

Version : 3 - 26/05/2026