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
INFORMATION FABRICANT / MANUFACTURER INFORMATION

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20, avenue du Stade de France - 93218 LA PLAINE SAINT-DENIS Cedex - Tel. (+33) 6 62 99 46 76 - efs.reactifs.hotline@efs.sante.fr

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For technical assistance or more information, the EFS Reactifs hotline is opened Monday to Friday from 9:00 a.m. to 5:00 p.m

Site de production / Manufacturing site

EFS Hauts de France – Normandie

10/12 Boulevard de Belfort - F-59000 LILLE cedex - Tél : (+33) 3 28 54 22 67 ou (+33) 3 28 54 49 31 - hfnio.upr-lille@efs.sante.fr

1. INTENDED USE AND INDICATION

SeroQal ARCHI 7 is a ready-to-use, unassigned quality control for the markers HBs Ag, HbC, Ab, HCV Ab, HIV-1 Ab, HTLV-I Ab, Syphilis and CMV Ab. It is available without an assigned value and is ready to use.

It is intended for in vitro use with Abbott ARCHITECT reagents and analyzers for serological testing in human plasma / serum. With its low reactivity, it allows for qualitative monitoring of analytical process performance and the detection of procedural errors or potential failures.

SeroQal ARCHI 7 is an independent control. It must not be substituted for the mandatory control reagents provided with the manufactured test kits.

SeroQal ARCHI 7 is an in vitro diagnostic medical device (IVDMD) in accordance with EU Regulation 2017/746.

2. COMPOSITION ET PACKAGING

Composition	Description	Packaging	Number of tests
Human plasma and serum collected from blood donors	Positive control for HBs Ag, HbC Ab, HCV Ab, HIV-1 Ab, HTLV-I Ab, Syphilis and CMV Ab markers	1 tube 3mL/tube	The sample volume allows at least 1 test for each marker

The human sera and plasma composing SeroQal are controlled with respect to markers HBs Ag, HbC Ab, HCV Ab, HIV p24 Ag, HIV-1 and HIV-2 Ab, HTLV-I and HTLV-II Ab and Syphilis Ab. Samples showing positive serology are treated with methylene blue to inactivate pathogens and thus reduce the risk of infection.

SeroQal are filtered to a 0.2µm (clarifying filtration) and are free of preservatives and antibiotics. They are non sterile.

3. DELIVERY AND STORAGE CONDITIONS

SeroQAL must be delivered at a temperature $\leq -20^{\circ}\text{C}$ (controlled temperature or dry ice transport).

Before opening, SeroQal must be placed and stored upon receipt in an enclosure $\leq -20^{\circ}\text{C}$. They are then stable until the expiration date indicated on the label.

After opening, SeroQal can be stored between $+2^{\circ}\text{C}$ and $+8^{\circ}\text{C}$ for 7 days or can be refrozen after use without exceeding 3 freeze-thaw cycles. Monitoring the number of cycles is the responsibility of the user.

4. PRECAUTIONS AND LIMITS FOR USE

SeroQal ARCHI 7 use is strictly reserved for health professionals and more specifically, laboratory technicians trained and authorized to perform analyses on the ARCHITECT analysers and to handle of biological samples.

Bad transport, storage or use conditions can lead erroneous results. It is recommended not to use and dispose of products if conditions are not guaranteed.

Do not use if tubes are damaged or if products show physical alteration indicating possible instability or degradation of the product.

5. SAFETY INSTRUCTIONS

SeroQal do not contain a hazardous ingredient, no substances with European Community workplace exposure limit values, or substances of very high concern (SVHC) in quantities greater than their respective reporting limit values. Consequently, the product is not supplied with a safety data sheet because the concentrations of chemicals contained in this product do not exceed the thresholds defined in the European Regulation 1907/2006.

Products are composed of blood products of human origin from donors subjected to mandatory blood screening tests according to French regulation. The wearing of personal protective equipment is recommended when handling the product.

SeroQal are not subject to the regulations of Transport of Dangerous Goods with regard to the ADR and IATA except for UN1845 labelling for shipment in dry ice.

SeroQal as well as all media that have been in contact with SeroQal must be handled as potentially infectious products. Disposal conditions and disinfection procedures must comply with current regulations.

6. PRINCIPLE OF USE

SeroQal ARCHI 7 was designed for use in each series of analyses under the same conditions as those applied to the biological samples.

It is formulated to produce a positive reactivity, the signal strength of which may vary depending on the laboratory, analyzer, and reagent lot used. It has no assigned value.

7. MATERIALS, REAGENT AND CONSUMABLES

The performance of the SeroQal ARCHI 7 is guaranteed with the use of the ARCHITECT analytical systems in combination with the materials, reagents and consumables specified in this IFU or validated by the user. The use of a different combination must be validated by the user.

Architect HBsAg - Anti-HBc - Anti-HCV - HIV Ag/Ab Combo - HTLV-I/II - Syphilis - CMV

8. PROTOCOL

The samples are ready for use. They must be left at room temperature until completely thawed, then **gently homogenized, avoiding foaming**.

The manufacturer's Abbott recommendations for biological samples to be analyzed tests procedure apply to SeroQal ARCHI 7.

The tube is labeled with a bar code for computerized lot number reading.

9. RESULTS INTERPRETATION

A POSITIVE result is expected:

- $0,8 \leq S/Co \leq 5,0$ for HBsAg, HbC Ab, HCV Ab, HIV-1 Ab, HTLV-I Ab and SYPHILIS Ab
- Concentration ≥ 6 AU/mL for CMV Ab.

If the results are not in accordance with the expected results, it may indicate unsatisfactory performance of the analytical process. It is then advised to implement the non-conformities management procedure of the user laboratory, to analyze the possible causes of this result and to invalidate the series if necessary.

10. PERFORMANCES

Performance studies were conducted on ARCHITECT i2000 analyzers. Reproducibility and repeatability analyses were performed at 4 sites using 3 different lots of SeroQal ARCHI 7. Repeatability: The 3 lots were tested in replicates of 10 times in the same series of analyses. The results were in accordance with the specifications announced for each of the serological markers. Reproducibility: The 3 lots were tested in replicates of 2 to 4 times per day for 2 to 7 days depending on the sites. All values ($n = 50$ to 71 depending on the serological markers studied) were in accordance with the specifications announced for each of the serological markers.

11. DETECTION OF ANOMALIES

Any serious incident occurring in connection with the device must be notified to the manufacturer and to the competent authority of the Member State in which the user is established.

Any incident related to the device's performance and/or safety must be reported to the manufacturer as soon as possible.