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1.- INTENDED USE AND INDICATIONS

NatQal B19-L is a ready to use Human Parvovirus B19 (B19V) DNA - positive quality control sample without any assigned value.

It is intended to be used with the cobas® DPX assay and the cobas® 6800/8800 systems manufactured by Roche Diagnostics, or with the Procleix® Parvo HAV Assay and the Tigris® and Panther® systems manufactured by Grifols, for Nucleic Acid Testing (NAT) in human plasma/serum.

NatQal B19-L is an in vitro medical device (IVMD) in accordance with EU Regulation 2017/746. This control sample is used to monitor performance of routine NAT analytical process. It allows the user to detect procedural errors or potential failures in the analytical process.

This device must be used under the same conditions as those applied to biological samples. It must not be substituted for the mandatory control reagents provided with the manufactured test kits.

2.- COMPOSITION AND PACKAGING

NatQal B19-L is manufactured from:

- ▶ human plasma positive for B19V DNA sourced from a blood donation, non-reactive for HIV-1, HCV and HAV RNA, HBV DNA, HBs antigen, anti-HIV-1/2, anti-HCV and anti-HTLV-1/2 antibodies;
- ▶ pooled human plasma sourced from blood donations collected, qualified and prepared in accordance with the good transfusion practices¹, filtered and tested non-reactive for HIV-1, HCV and HAV RNA and HBV and B19V DNA.

NatQal B19-L is supplied in two different packaging sizes:

- Reference **247901**: NatQal B19-L - 1 ml x 50 vials packaging

- Reference **247900**: NatQal B19-L - 1 ml x 10 vials packaging

The product barcode, which allows computer-based reading of the batch code, is included inside the packaging.

3.- DELIVERY AND STORAGE CONDITIONS

NatQal B19-L is carried at controlled temperature allowing the maintenance of a temperature $\leq -20^{\circ}\text{C}$. Store NatQal B19-L upon receipt at -20°C or colder until use or expiration date.

To prevent leakage, store vials of NatQal B19-L upright.

INDICATIONS OF REAGENT INSTABILITY OR DETERIORATION

NatQal B19-L is delivered and stored in frozen form. After thawing and homogenization, it appears as a translucent yellow liquid.

Alterations in physical appearance may indicate instability or deterioration of NatQal B19-L. In this case, do not use the device and please contact EFS Réactifs hotline for assistance.

4.- WARNINGS AND PRECAUTIONS – LIMITATIONS**CAUTION**

Handle NatQal B19-L and all human blood derivatives **as though capable of transmitting infectious agents** according to Good Laboratory Practices.

SAFETY PRECAUTIONS

Only health professionals, and especially laboratory technicians according to the French Public Health Code L4352-1 to 9 or local legislation, qualified to perform qualitative and/or quantitative NAT assays should use this device. It must be handled in appropriate containment facilities by competent personnel who are fully trained in the handling of infectious products and in accordance with national and local safety rules.

Rules for the handling of contaminated human blood products must be applied². Use the World Health Organization (WHO) recommended universal precautions for handling NatQal B19-L and human blood.

Do not eat or drink in areas where specimens are being handled. Do not pipette by mouth. Wear appropriate laboratory protections (lab coats, gloves, goggles, etc.).

Clean any spillage by immediately wiping up with 0.5% sodium hypochlorite solution.

Any serious incident related to the device must be notified to the manufacturer and to the competent authority of the Member State in which the user is established.

After use, dispose of all specimens, controls and materials used in testing as though they contain infectious agents, in accordance with laboratory procedures and Good Laboratory Practices.

HANDLING PRECAUTIONS

Do not use NatQal B19-L beyond the expiration date.

Avoid contaminations when opening and closing vials.

Do not transfer the product to another recipient.

LIMITS OF USE

NatQal B19-L must be stored and used in accordance with the instructions given in this document. Test procedures and interpretation of results provided by test kits manufacturers must be followed. Deviations from procedures recommended by test kits manufacturers may produce unreliable results.

The performance characteristics of NatQal B19-L have been established only for the detection of B19V DNA.

Performances of NatQal B19-L are not guaranteed in case of inappropriate shipping, storage conditions or use of expired products.

INFORMATION FABRICANT / MANUFACTURER INFORMATION**Fabricant légal / Manufacturer**

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For technical assistance, information or report, EFS Réactifs hotline is open Monday to Friday from 9:00 am to 5:00 pm.

5.- USE PRINCIPLE

Routine use of independent positive controls enables laboratories to monitor day-to-day and operator variation, lot-to-lot performance of test kits, and can assist the user in identifying errors and failures, and in revealing possible defects in the analytical chain or in the reagents used for the detection of B19V DNA.

NatQal B19-L simulates a positive sample or a positive sample pool. It follows all the steps of the analytical process going from the preparation of the sample, to its amplification and detection.

In this context, NatQal B19-L participates in the validation of a given series of assays, as well as for the validation upon receipt of the reagents used.

6.- MATERIAL, REAGENTS AND OTHER CONSUMABLES

NatQal B19-L has been validated on analytical systems supplied by Grifols (Procleix® Tigris® system and Procleix® Panther® system – Procleix® Parvo HAV Assay) and by Roche Diagnostics (cobas® 6800/8800 System – cobas® DPX Assay).

Any replacement of a reagent/material combination must be subjected to method validation by the user.

7.- PROCEDURE**MATERIALS REQUIRED BUT NOT PROVIDED**

Refer to user manuals provided by the manufacturers of the NAT kits used.

INSTRUCTIONS FOR USE

NatQal B19-L **must be quickly and completely thawed at 37°C** to avoid the formation of cryoprecipitates. Gently mix the tube at this stage until ice clot has disappeared.

Briefly vortex the tube before use. If necessary, remove residual drops of solution in the stopper and open the vial carefully. If necessary, remove bubbles or foam using a sterile pipette.

Do not centrifuge the tube.

NatQal B19-L is for single use, and should not be reused.

After thawing, the vial should be used immediately.

NatQal B19-L should be included in any test run in accordance with the procedure provided by the manufacturer for unknown samples. It must not be substituted for the control reagents provided with manufactured test kits.

This device is ready-to-use. It has been designed to be integrated directly into the analytical chain at the beginning of the sample preparation step or after the pooling step if it is processed. It does not undergo the pooling step.

8.- EXPECTED RESULTS

NatQal B19-L should yield a **POSITIVE result**, in accordance with experimental data and the manufacturers' instructions for use of NAT reagents for B19V DNA screening, with a viral load targeted at 5,000 IU/mL, i.e., $3.7 \pm 0.3 \log_{10}$ IU/mL.

If the results are not consistent with those expected, this may indicate an unsatisfactory performance of the analytical process. The user is then advised to proceed to the analysis of the possible causes, to invalidate the entire run and to inform EFS Réactifs manufacturer if necessary.

9.- DIAGNOSTIC PERFORMANCES

Analytical performances and stability studies were carried out on the analytical systems of the company Grifols (Procleix® Tigris® system and Procleix® Panther® system – Procleix® Parvo HAV Assay) as well on the analytical systems of the company Roche Diagnostics (cobas® 6800/8800 System – cobas® DPX Assay). The results obtained are consistent with data from NAT reagents' manufacturers for the detection of B19V DNA. 100% of NatQal B19-L have been tested POSITIVE.

The viral load of the sample is targeted at 5,000 IU/mL for B19V DNA, i.e., $3.7 \log_{10}$ IU/mL. For information purposes, the individual values obtained during performance studies ranged from 3.5 to 4.5 \log_{10} IU/mL ($\text{CI}_{95\%}$).

Procedures for implementing a quality assurance program and monitoring test performance on a routine basis must be established by each individual laboratory.

10.- MATERIAL SAFETY SHEET

NatQal B19-L does not contain any hazardous chemical products, nor "substances of concern", nor "substances of very high concern" (SVHC) requiring limit values for exposure monitoring for the workplace and/or the environment according to European Community. Therefore, no safety data sheet is provided because the chemicals concentrations contained in this product do not exceed the thresholds defined in the European Regulation 1907/2006.

11.- LATEST UPDATES

Version / Date	Update
Ve 1 / 22/05/2026	Alignment with IVDR requirements Minor text corrections Layout review Details about the intended use of the device Modification of the manufacturer information box

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