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
1.- INTENDED USE AND INDICATIONS

NatQal HEV is a ready to use Hepatitis E virus (HEV) RNA - positive quality control sample without any assigned value.

It is intended to be used with the cobas® HEV assay and the cobas® 6800/8800 systems manufactured by Roche Diagnostics, or with the Procleix® HEV Assay or Procleix® UltrioPlex E Assay and the Panther® systems manufactured by Grifols, for Nucleic Acid Testing (NAT) in human plasma/serum.

NatQal HEV is an in vitro medical device (IVMD) in accordance with EU Regulation 2017/746. This control sample is used to monitor performance of routine NAT analytical process. It allows the user to detect procedural errors or potential failures in the analytical process.

This device must be used under the same conditions as those applied to biological samples. It must not be substituted for the mandatory control reagents provided with the manufactured test kits.

2.- COMPOSITION AND PACKAGING

NatQal HEV is manufactured from:

- human plasma sourced from a 3c genotype HEV RNA blood donation, non-reactive for HIV-1 and HCV RNA, HBV DNA, HBs antigen, anti-HIV-1/2, anti-HCV and anti-HTLV-1/2 antibodies;
- pooled human plasma sourced from blood donations collected, qualified and prepared in accordance with the good transfusion practices¹, filtered and tested non-reactive for HIV-1, HCV and HEV RNA and HBV DNA.

NatQal HEV is supplied in two different packaging sizes:

- Reference **274284**: NatQal HEV - 1 ml x 50 vials packaging
- Reference **274283**: NatQal HEV - 1 ml x 10 vials packaging

The product barcode, which allows computer-based reading of the batch code, is included inside the packaging.

3.- DELIVERY AND STORAGE CONDITIONS

NatQal HEV is carried at controlled temperature allowing the maintenance of a temperature ≤ -20°C. Store NatQal HEV upon receipt at -20°C or colder until use or expiration date.

To prevent leakage, store vials of NatQal HEV upright.

INDICATIONS OF REAGENT INSTABILITY OR DETERIORATION

NatQal HEV is delivered and stored in frozen form. After thawing and homogenization, it appears as a translucent yellow liquid.

Alterations in physical appearance may indicate instability or deterioration of NatQal HEV. In this case, do not use the device and please contact EFS Réactifs hotline for assistance.

4.- WARNINGS AND PRECAUTIONS – LIMITATIONS**CAUTION**

Handle NatQal HEV and all human blood derivatives **as though capable of transmitting infectious agents** according to Good Laboratory Practices.

SAFETY PRECAUTIONS

Only health professionals, and especially laboratory technicians according to the French Public Health Code L4352-1 to 9 or local legislation, qualified to perform qualitative and/or quantitative NAT assays should use this device. It must be handled in appropriate containment facilities by competent personnel who are fully trained in the handling of infectious products and in accordance with national and local safety rules.

Rules for the handling of contaminated human blood products must be applied². Use the World Health Organization (WHO) recommended universal precautions for handling NatQal HEV and human blood.

Do not eat or drink in areas where specimens are being handled. Do not pipette by mouth. Wear appropriate laboratory protections (lab coats, gloves, goggles, etc.).

Clean any spillage by immediately wiping up with 0.5% sodium hypochlorite solution.

Any serious incident related to the device must be notified to the manufacturer and to the competent authority of the Member State in which the user is established.

After use, dispose of all specimens, controls and materials used in testing as though they contain infectious agents, in accordance with laboratory procedures and Good Laboratory Practices.

HANDLING PRECAUTIONS

Do not use NatQal HEV beyond the expiration date.

Avoid contaminations when opening and closing vials.

Do not transfer the product to another recipient.

LIMITS OF USE

NatQal HEV must be stored and used in accordance with the instructions given in this document. Test procedures and interpretation of results provided by test kits manufacturers must be followed. Deviations from procedures recommended by test kits manufacturers may produce unreliable results.

The performance characteristics of NatQal HEV have been established only for the detection of HEV RNA.

Performances of NatQal HEV are not guaranteed in case of inappropriate shipping, storage conditions or use of expired products.

INFORMATION FABRICANT / MANUFACTURER INFORMATION**Fabricant légal / Manufacturer**

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For technical assistance, information or report, EFS Réactifs hotline is open Monday to Friday from 9:00 am to 5:00 pm.

5.- USE PRINCIPLE

Routine use of independent positive controls enables laboratories to monitor day-to-day and operator variation, lot-to-lot performance of test kits, and can assist the user in identifying errors and failures, and in revealing possible defects in the analytical chain or in the reagents used for the detection of HEV RNA.

NatQal HEV simulates a positive sample or a positive sample pool. It follows all the steps of the analytical process going from the preparation of the sample, to its amplification and detection.

In this context, NatQal HEV participates in the validation of a given series of assays, as well as for the validation upon receipt of the reagents used.

6.- MATERIAL, REAGENTS AND OTHER CONSUMABLES

NatQal HEV has been validated on analytical systems supplied by Grifols (Procleix® Panther® System - Procleix® HEV Assay or Procleix® UltrioPlex E Assay) or Roche Diagnostics (cobas® 6800/8800 – cobas® HEV Assay).

Any replacement of a reagent/material combination must be subjected to method validation by the user.

7.- PROCEDURE**MATERIALS REQUIRED BUT NOT PROVIDED**

Refer to user manuals provided by the manufacturers of the NAT kits used.

INSTRUCTIONS FOR USE

NatQal HEV **must be quickly and completely thawed at 37°C** to avoid the formation of cryoprecipitates. Gently mix the tube at this stage until ice clot has disappeared.

Briefly vortex the tube before use. If necessary, remove residual drops of solution in the stopper and open the vial carefully. If necessary, remove bubbles or foam using a sterile pipette.

Do not centrifuge the tube.

NatQal HEV is for single use, and should not be reused.

After thawing, the vial should be used immediately.

NatQal HEV should be included in any test run in accordance with the procedure provided by the manufacturer for unknown samples. It must not be substituted for the control reagents provided with manufactured test kits.

This device is ready-to-use. It has been designed to be integrated directly into the analytical chain at the beginning of the sample preparation step or after the pooling step if it is processed. It does not undergo the pooling step.

8.- EXPECTED RESULTS

NatQal HEV should give a **POSITIVE result**, in accordance with experimental data and the manufacturers' instructions for use of NAT reagents for HEV RNA screening, for the interpretation of results for a reactive sample.

If the results are not consistent with those expected, this may indicate an unsatisfactory performance of the analytical process. The user is then advised to proceed to the analysis of the possible causes, to invalidate the entire run and to inform EFS Réactifs manufacturer if necessary.

9.- DIAGNOSTIC PERFORMANCES

Analytical performances and stability studies were carried out on the analytical systems of the company Grifols (Procleix® Panther® System - Procleix® HEV Assay or Procleix® UltrioPlex E Assay) as well on the analytical systems of the company Roche Diagnostics (cobas® 6800/8800 System – cobas® HEV Assay). The results obtained are consistent with data from NAT reagents' manufacturers for the detection of HEV RNA. 100% of NatQal HEV have been tested POSITIVE.

The viral load of the sample is estimated as an indication at 60 IU/ml. This viral load represents at least thrice the 95% detection threshold (SD95%) of detection techniques³, determined routinely by the French Blood Institute (EFS) laboratories.

Procedures for implementing a quality assurance program and monitoring test performance on a routine basis must be established by each individual laboratory.

10.- MATERIAL SAFETY SHEET

NatQal HEV does not contain any hazardous chemical products, nor "substances of concern", nor "substances of very high concern" (SVHC) requiring limit values for exposure monitoring for the workplace and/or the environment according to European Community. Therefore, no safety data sheet is provided because the chemicals concentrations contained in this product do not exceed the thresholds defined in the European Regulation 1907/2006.

11.- LATEST UPDATES

Version / Date	Update
Ve 3 / 16/01/2025	Replacement of trade name RC DGV MONO HEV by NatQal HEV Update of the hotline contact information Replacement of UPR by EFS Réactifs
Ve 4 / 22/12/2025	Minor text corrections Layout review Modification of the manufacturer information box Modification of the product reference table in the footer

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